

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-PBRM1



Numéro de catalogue: 12563-1-AP

2 Publications

Informations de base

Numéro de catalogue:

12563-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop and 280 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG3148

Numéro d'acquisition GenBank:

BC015323

Identification du gène (NCBI):

55193

Nom complet:

polybromo 1

MW calculé

180 kDa

MW observés:

180-195 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Demandes citées:

Dot blot, IP, WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain

Contrôles positifs:

WB: cellules Jurkat,

Informations générales

PBRM1, also named hPB1 BAF180 and Polybromo-1D, acts as a negative regulator of cell proliferation. PBRM1 is the second major mutated clear cell RCC cancer gene after VHL. PBRM1 loss is associated with low grade features.

Publications notables

Autrice	Pubmed ID	Journal	Application
Sarah E Hopson	30009957	Protein Expr Purif	Dot blot
Wentao Liu	35368029	Cell Death Dis	WB, IP

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Jurkat cells were subjected to SDS PAGE followed by western blot with 12563-1-AP (PBRM1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.