

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-IL-18BP



Numéro de catalogue: 14153-1-AP

## Informations de base

Numéro de catalogue:

14153-1-AP

Taille:

150ul, Concentration: 300 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG5335

Numéro d'acquisition GenBank:

BC044215

Identification du gène (NCBI):

10068

Nom complet:

interleukin 18 binding protein

MW calculé

194 aa, 21 kDa

MW observés:

21 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

## Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, rat, souris

Contrôles positifs:

WB : cellules A549, cellules HeLa

## Informations générales

IL18BP, also named as adekinig-alfa, has 3 isoforms (Isoform A 21kd; B 12kd and D 17kd). Isoform A binds to IL-18 and inhibits its activity. IL18BP functions as an inhibitor of the early TH1 cytokine response. Highly glycosylated of IL18BP will migrate to be 40kd. A 58kd protein is IL18BP crosslinked IL18.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 14153-1-AP (IL-18BP antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.