

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-SERPINB9



Numéro de catalogue: 14910-1-AP

Informations de base

Numéro de catalogue:

14910-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop and 307 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG6707

Numéro d'acquisition GenBank:

BC002538

Identification du gène (NCBI):

5272

Nom complet:

serpin peptidase inhibitor, clade B (ovalbumin), member 9

MW calculé

42 kDa

MW observés:

42 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:1000-1:3000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules K-562, tissu placentaire humain

Informations générales

SERPINB9, also named as PI-9, belongs to the serpin family. SERPINB9 is an intracellular inhibitor of granzyme B (grB) that protects cytotoxic lymphocytes from grB-mediated death. It has a nucleo-cytoplasmic distribution and is produced in CD8+ T cells and natural killer cells (PMID: 28024184). Also, SERPINB9 is expressed by vascular SMCs and endothelial cells to avoid GZB-induced apoptosis (PMID: 28743062). The molecular weight of SERPINB9 is 42 kDa.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 14910-1-AP (SERPINB9 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.