

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-PDCL



Numéro de catalogue: 16057-1-AP

Informations de base

Numéro de catalogue:	16057-1-AP	Numéro d'acquisition GenBank:	BC017227	Méthode de purification:	Purification par affinité contre l'antigène
Taille:	150ul, Concentration: 400 µg/ml by Nanodrop;	Identification du gène (NCBI):	5082	Dilutions recommandées:	WB 1:1000-1:6000 IHC 1:50-1:500
Hôte:	Lapin	Nom complet:	phosducin-like		
Isotype:	IgG	MW calculé	301 aa, 34 kDa		
Immunogen Catalog Number:	AG9029	MW observés:	40-42 kDa		

Applications

Applications testées:
IHC, WB, ELISA
Spécificité de l'espèce:
Humain, rat, souris

Contrôles positifs:

WB : cellules HEK-293, cellules Jurkat, cellules NCI-H1299

IHC : tissu placentaire humain, tissu cérébral de souris

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9.0; (*) À défaut, 'le démasquage de l'antigène peut être 'effectué avec un tampon citrate pH 6,0.

Informations générales

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




Various lysates were subjected to SDS PAGE followed by western blot with 16057-1-AP (PDCL antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human placenta using 16057-1-AP (PDCL antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human placenta using 16057-1-AP (PDCL antibody) at dilution of 1:50 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 16057-1-AP (PDCL antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).