

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-MAGEC1-Specific



Numéro de catalogue: 16176-1-AP

## Informations de base

Numéro de catalogue:

16176-1-AP

Taille:

150ul, Concentration: 473 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_005462

Identification du gène (NCBI):

9947

Nom complet:

melanoma antigen family C, 1

MW calculé

124 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain

## Informations générales

This gene is a member of the melanoma antigen gene (MAGE) family. The proteins of this family are tumor-specific antigens. Expression of this protein was not seen in normal tissues, with the exception of testis. This protein contains a large number of unique short repetitive sequences in front of the MAGE-homologous sequence, and therefore is about 800 aa longer than the other MAGE proteins. (From OMIM)

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

