

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-KDEL<sup>C</sup>1



Numéro de catalogue:17018-1-AP

## Informations de base

Numéro de catalogue:

17018-1-AP

Taille:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG9999

Numéro d'acquisition GenBank:

BC001297

Identification du gène (NCBI):

79070

Nom complet:

KDEL (Lys-Asp-Glu-Leu) containing 1

MW calculé

58 kDa

MW observés:

50 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

WB,ELISA

Spécificité de l'espèce:

Humain

## Informations générales

KDEL<sup>C</sup>1, also named as EP58 and ERp58, belongs to the KDEL<sup>C</sup> family. Soluble proteins in the endoplasmic reticulum (ER) contain a specific carboxyterminal sequence KDEL (Lys-Asp-Glu-Leu) and include the coat proteins required for vesicle budding from the ER. KDEL<sup>C</sup>1 shares no significant homology to any of the known ER-resident proteins. It is primarily expressed in embryo, placenta and adult heart.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 17018-1-AP (KDEL1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.