

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-XAF1



Numéro de catalogue: 17837-1-AP

Informations de base

Numéro de catalogue:

17837-1-AP

Taille:

150ul

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG12110

Numéro d'acquisition GenBank:

BC073156

Identification du gène (NCBI):

54739

Nom complet:

XIAP associated factor 1

MW calculé

301 aa, 35 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

IF, ELISA

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of HeLa cells using 17837-1-AP (XAF1 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).