

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-LAMC2



Numéro de catalogue: 19698-1-AP

3 Publications

## Informations de base

Numéro de catalogue:

19698-1-AP

Taille:

150ul, Concentration: 900 µg/ml by Nanodrop and 413 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_005562

Identification du gène (NCBI):

3918

Nom complet:

laminin, gamma 2

MW calculé

131 kDa

MW observés:

140 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

## Applications

Applications testées:

WB, ELISA

Demandes citées:

IHC, WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain, rat

Contrôles positifs:

WB : cellules A431,

## Informations générales

LAMC2, also named as LAMB2T and LAMNB2, is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components. LAMC2 binds to cells via a high affinity receptor. Laminin exerts cell-scattering activity toward a wide variety of cells, including epithelial, endothelial, and fibroblastic cells.

## Publications notables

| Autrice         | Pubmed ID | Journal        | Application |
|-----------------|-----------|----------------|-------------|
| Giricea Lorusso | 36070364  | Sci Transl Med | WB, IHC     |
| Caiqun Huang    | 33179081  | Mol Med Rep    | WB          |
| Caiming Xu      | 32244009  | J Proteomics   | WB          |

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



A431 cells were subjected to SDS PAGE followed by western blot with 19698-1-AP (LAMC2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.