

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-HTR1B



Numéro de catalogue: 19960-1-AP

## Informations de base

Numéro de catalogue:

19960-1-AP

Taille:

150ul

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_000863

Identification du gène (NCBI):

3351

Nom complet:

5-hydroxytryptamine (serotonin)

receptor 1B

MW calculé

44 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

IHC, WB, ELISA

Spécificité de l'espèce:

Humain, rat, souris

## Informations générales

HTR1B, also named as 5-HT1B;5-HT1DB;HTR1D2;HTR1DB and 512, belongs to the G-protein coupled receptor 1 family. HTR1B is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of HTR1B is mediated by G proteins that inhibit adenylate cyclase activity.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

