

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-TRERF1



Numéro de catalogue: 20062-1-AP

1 Publications

Informations de base

Numéro de catalogue:

20062-1-AP

Taille:

150ul, Concentration: 450 µg/ml by Nanodrop and 327 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM_033502

Identification du gène (NCBI):

55809

Nom complet:

transcriptional regulating factor 1

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:200-1:1000

MW calculé

132 kDa

MW observés:

132-140 kDa

Applications

Applications testées:

WB, ELISA

Demandes citées:

IHC

Spécificité de l'espèce:

Humain, souris

Espèces citées:

Humain

Contrôles positifs:

WB : cellules MCF7, cellules MCF-7

Informations générales

TRERF1, also named as BCAR2, RAPA and TREP132, activates transcription of CYP11A1. TRERF1 interaction with CREBBP and EP300 results in a synergistic transcriptional activation of CYP11A1. The antibody recognizes the C-term of TRERF1.

Publications notables

Autrice	Pubmed ID	Journal	Application
Jianlong Wang	35733512	Cancer Manag Res	IHC

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



MCF7 cells were subjected to SDS PAGE followed by western blot with 20062-1-AP (TRERF1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.