

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-HDAC9



Numéro de catalogue: 23328-1-AP

Informations de base

Numéro de catalogue:

23328-1-AP

Numéro d'acquisition GenBank:

BC152405

Méthode de purification:

Purifié par affinité contre l'antigène

Taille:

150ul, Concentration: 433 µg/ml by Bradford method using BSA as the standard;

Identification du gène (NCBI):

9734

Nom complet:

histone deacetylase 9

Hôte:

Lapin

MW calculé

1011 aa, 111 kDa

Isotype:

IgG

MW observés:

100 kDa

Immunogen Catalog Number:

AG19930

Applications

Applications testées:

IHC, ELISA

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 23328-1-AP (HDAC9 Antibody) at dilution of 1:200 (under 40x lens).