

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNF253



Numéro de catalogue: 25152-1-AP

Informations de base

Numéro de catalogue:

25152-1-AP

Taille:

150ul, Concentration: 1300 µg/ml by Nanodrop and 627 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG18570

Numéro d'acquisition GenBank:

BC125063

Identification du gène (NCBI):

56242

Nom complet:

zinc finger protein 253

MW calculé

499 aa, 58 kDa

MW observés:

58-65 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:200-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules HepG2, cellules SMMC-7721

Informations générales

ZNF253, also named as Zinc finger protein 411 or BMZF1, is a 499 amino acid protein, which contains one KRAB domain and eleven C2H2-type zinc fingers. ZNF253 belongs to the krueppel C2H2-type zinc-finger protein family and is expressed in bone marrow and in monocytic and immature erythroid cell lines. ZNF253 may function as a transcription factor.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HepG2 cells were subjected to SDS PAGE followed by western blot with 25152-1-AP (ZNF253 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.