

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNF711



Numéro de catalogue: 25166-1-AP

Informations de base

Numéro de catalogue:

25166-1-AP

Taille:

150ul, Concentration: 650 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG18309

Numéro d'acquisition GenBank:

BC067294

Identification du gène (NCBI):

7552

Nom complet:

zinc finger protein 711

MW calculé

761 aa, 86 kDa

MW observés:

115 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, souris

Contrôles positifs:

WB : cellules SKOV-3, cellules Jurkat, tissu testiculaire de souris

Informations générales

ZNF711 also named as CMPX1 or Zinc finger protein 6 is a 761 amino acid protein, which contains 12 C2H2-type zinc fingers and belongs to the krueppel C2H2-type zinc finger protein family. ZNF711 localizes in the nucleus and is expressed in neural tissue. ZNF711 as a transcription factor is required for brain development. ZNF711 mutation may have an association with mental retardation.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



SKOV-3 cells were subjected to SDS PAGE followed by western blot with 25166-1-AP (ZNF711 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



mouse testis tissue were subjected to SDS PAGE followed by western blot with 25166-1-AP (ZNF711 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.