

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNF639



Numéro de catalogue: 25817-1-AP

Informations de base

Numéro de catalogue:
25817-1-AP

Taille:
150ul, Concentration: 300 µg/ml by
Nanodrop and 200 µg/ml by Bradford
method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG19854

Numéro d'acquisition GenBank:
BC026181

Identification du gène (NCBI):
51193

Nom complet:
zinc finger protein 639

MW calculé
485 aa, 56 kDa

MW observés:
65-70 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:200-1:1000
IP 0.5-4.0 ug for IP and 1:200-1:1000
for WB

Applications

Applications testées:
IP, WB, ELISA

Spécificité de l'espèce:
Humain

Contrôles positifs:

WB : cellules HeLa,

IP : cellules HeLa,

Informations générales

ZNF639, also named as ZASC1, is a 485 amino acid protein, which locates in the nucleus and contains 8 C2H2-type zinc fingers and belongs to the krueppel C2H2-type zinc-finger protein family. ZNF639 binds DNA and may function as a transcriptional repressor.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HeLa cells were subjected to SDS PAGE followed by western blot with 25817-1-AP (ZNF639 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.

IP Result of anti-ZNF639 (IP:25817-1-AP, 4ug; Detection:25817-1-AP 1:300) with HeLa cells lysate 3960ug.