

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-POU5F2



Numéro de catalogue: 26362-1-AP

Informations de base

Numéro de catalogue:

26362-1-AP

Taille:

150ul, Concentration: 650 µg/ml by Nanodrop and 367 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG23803

Numéro d'acquisition GenBank:

BC029532

Identification du gène (NCBI):

134187

Nom complet:

POU domain class 5, transcription factor 2

MW observés:

42-45 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:200-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB: cellule NCCIT,

Informations générales

POU5F2, also named as POU domain, class 5, transcription factor 2m, is a 328 amino acid protein, which contains 1 homeobox DNA-binding domain and 1 POU-specific domain and belongs to the POU transcription factor family. POU5F2 as a transcription factor, localizes in the nucleus and may exert a regulatory function in meiotic events that are required for terminal differentiation of male germ cell.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



NCCIT cell were subjected to SDS PAGE followed by western blot with 26362-1-AP (POU5F2 Antibody) at dilution of 1:200 incubated at room temperature for 1.5 hours.