

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZBTB7B



Numéro de catalogue: 28127-1-AP

Informations de base

Numéro de catalogue: 28127-1-AP	Numéro d'acquisition GenBank: BC012070	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul, Concentration: 600 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 51043	Dilutions recommandées: WB 1:500-1:2000
Hôte: Lapin	Nom complet: zinc finger and BTB domain containing 7B	
Isotype: IgG	MW calculé 539 aa, 58 kDa	
Immunogen Catalog Number: AG27970	MW observés: 70 kda	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules HepG2, cellules HeLa
Spécificité de l'espèce: Humain, souris	

Informations générales

ZBTB7B, also known as Zfp-67 or Zinc finger protein Th-POK, is a 539 amino acid protein, which localizes in nucleus. ZBTB7B as a transcription regulator that acts as a key regulator of lineage commitment of immature T-cell precursors. ZBTB7B is necessary and sufficient for commitment of CD4 lineage, while its absence causes CD8 commitment. The calculated molecular weight of ZBTB7B is about 58 kDa, but modified protein is 70 kDa.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 28127-1-AP (ZBTB7B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.