

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-SREBF2



Numéro de catalogue: 28212-1-AP

Phare

34 Publications

Informations de base

Numéro de catalogue:

28212-1-AP

Taille:

150ul, Concentration: 900 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG28205

Numéro d'acquisition GenBank:

BC056158

Identification du gène (NCBI):

6721

Nom complet:

sterol regulatory element binding transcription factor 2

MW calculé

124 kDa

MW observés:

73-75 kDa, 120-130 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:1000-1:8000

Applications

Applications testées:

WB, ELISA

Demandes citées:

IHC, WB

Spécificité de l'espèce:

Humain

Espèces citées:

Humain, souris, Hamster

Contrôles positifs:

WB : cellules HUVEC, cellules HeLa, cellules HepG2

Informations générales

SREBF2, also named as BHLHD2, is a 1141 amino acid protein, which contains 1 bHLH domain and belongs to the SREBP family. SREBF2 localizes in the endoplasmic reticulum membrane and is ubiquitously expressed in adult and fetal tissues. SREBF2 as a transcriptional activator is required for lipid homeostasis. SREBF2 expression is dynamically regulated in response to nutritional and hormonal cues. SREBF2 exists two isoform with molecular weight 124 kDa and 73 kDa.

Publications notables

Autrice	Pubmed ID	Journal	Application
Shaofan Hu	36174386	Redox Biol	WB
Mingyue Tao	34518524	Cell Death Dis	WB
Hui-Yu Luo	36286043	Curr Issues Mol Biol	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 28212-1-AP (SREBF2 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.