

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-CCSAP



Numéro de catalogue: 28644-1-AP

Informations de base

Numéro de catalogue:

28644-1-AP

Taille:

150ul, Concentration: 220 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG23785

Numéro d'acquisition GenBank:

BC015419

Identification du gène (NCBI):

126731

Nom complet:

chromosome 1 open reading frame 96

MW observés:

30-36 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

rat, souris

Contrôles positifs:

WB : tissu cérébral de souris, tissu cérébral de rat

Informations générales

CCSAP (Centriole, cilia and spindle-associated protein) colocalizes with polyglutamylated tubulin to centrioles, spindle microtubules, and cilia in human tissue culture cells (PMID: 22493317). Plays a role in microtubule (MT) stabilization and this stabilization involves the maintenance of NUMA1 at the spindle poles. It can be detected as around 30-36 kDa.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 28644-1-AP (C1orf96 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.