

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-PPP2R5A



Numéro de catalogue: 29365-1-AP

Informations de base

Numéro de catalogue:	29365-1-AP	Numéro d'acquisition GenBank:	BC022474	Méthode de purification:	Purification par affinité contre l'antigène
Taille:	150ul, Concentration: 800 µg/ml by Nanodrop;	Identification du gène (NCBI):	5525	Dilutions recommandées:	WB 1:2000-1:10000
Hôte:	Lapin	Nom complet:	protein phosphatase 2, regulatory subunit B', alpha isoform		
Isotype:	IgG	MW calculé	486 aa, 56 kDa		
Immunogen Catalog Number:	AG31072	MW observés:	50-56 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : cellules Jurkat, tissu cardiaque de rat, tissu cardiaque de souris
Spécificité de l'espèce:	Humain, rat, souris		

Informations générales

Protein phosphatase 2, regulatory subunit B (B56), alpha isoform (PPP2R5A) is one of the regulatory subunits of the protein phosphatase 2A (PP2A), which is a major member of protein serine/threonine phosphatases in cells. PPP2R5A can regulate the cellular location, substrate specification and protein phosphatase function of PP2A, through which PPP2R5A plays an important role in many cellular activities. It has been reported that PPP2R5A relates with many diseases including cancers by regulating many crucial signaling pathways involved in P53, Bcl-2, CDK, MAPK, JAK/STAT, c-Myc and β -Catenin.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 29365-1-AP (PPP2R5A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.