

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-BRCA2



Numéro de catalogue: 29450-1-AP

3 Publications

Informations de base

Numéro de catalogue:

29450-1-AP

Taille:

150ul, Concentration: 1300 µg/ml by 675

Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG31064

Numéro d'acquisition GenBank:

NM_000059

Identification du gène (NCBI):

675

Nom complet:

breast cancer 2, early onset

MW calculé

384 kDa

MW observés:

384 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:1000-1:8000

Applications

Applications testées:

WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain

Espèces citées:

Humain

Contrôles positifs:

WB : cellules HeLa, cellules HepG2, cellules MCF-7

Informations générales

BRCA2, also named as FACD and FANCD1, is involved in double-strand break repair and/or homologous recombination. BRCA2 may participate in S phase checkpoint activation. Defects in BRCA2 are a cause of genetic susceptibility to breast cancer (BC) which also called susceptibility to familial breast-ovarian cancer type 2 (BROVCA2).

Publications notables

Autrice	Pubmed ID	Journal	Application
Taihao Chang	36576015	Prostate	WB
Nicole L Batenburg	37569794	Int J Mol Sci	WB
Yingwei Li	37391593	Oncogene	WB

Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 29450-1-AP (BRCA2 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.