

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-T7 RNA polymerase



Numéro de catalogue: 29943-1-AP

Informations de base

Numéro de catalogue:

29943-1-AP

Taille:

150ul, Concentration: 600 µg/ml by
Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG32001

Numéro d'acquisition GenBank:

NP_041960

Identification du gène (NCBI):

1261050

Nom complet:

T7 RNA polymerase

MW calculé

99 kDa

Méthode de purification:

Purification par affinité contre
l'antigène

Dilutions recommandées:

WB 1:2000-1:10000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Bacteriophage T7

Contrôles positifs:

WB : Protéine recombinante,

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Données de validation sélectionnées



Recombinant protein were subjected to SDS PAGE followed by western blot with 29943-1-AP (T7 RNA polymerase antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.