

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-DCBLD2



Numéro de catalogue: 30245-1-AP

Informations de base

Numéro de catalogue:	30245-1-AP	Numéro d'acquisition GenBank:	BC029658	Méthode de purification:	Purification par affinité contre l'antigène
Taille:	150ul, Concentration: 500 µg/ml by Nanodrop;	Identification du gène (NCBI):	131566	Dilutions recommandées:	WB 1:1000-1:6000
Hôte:	Lapin	Nom complet:	discoïdin, CUB and LCCL domain containing 2		
Isotype:	IgG	MW calculé	775 aa, 89 kDa		
Immunogen Catalog Number:	AG32823	MW observés:	120 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : cellules A549, cellules U2OS
Spécificité de l'espèce:	Humain		

Informations générales

DCBLD2, also known as ESDN or CLCP1, was first identified from human coronary arterial cells by a signal sequence trap method. DCBLD2 regulates proliferation, migration, and invasion in various tumors, including glioma, hypopharyngeal squamous cell carcinoma, and myxofibrosarcoma. The human DCBLD2 protein has a predicted molecular mass of 80 kDa based on its amino acid sequence, although it regularly displays an effective molecular mass of 130 kDa via SDS-PAGE (PMID: 30902898).

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 30245-1-AP (DCBLD2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.