

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-BLM



Numéro de catalogue: 30254-1-AP

Informations de base

Numéro de catalogue:

30254-1-AP

Taille:

150ul, Concentration: 550 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG32996

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

641

Nom complet:

Bloom syndrome, RecQ helicase-like

MW calculé

159 kDa

MW observés:

159 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:1000-1:4000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules NIH/3T3, cellules HEK-293T, cellules HeLa, cellules Jurkat, cellules K-562

Informations générales

BLM participates in DNA replication and repair, it is involved in 5'-end resection of DNA during double-strand break (DSB) repair. BLM can also be recruited by the KHDC3L-OOEP scaffold to DNA replication forks where it is retained by TRIM25 ubiquitination, it thereby promotes the restart of stalled replication forks.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 30254-1-AP (BLM antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.