

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-HLA-DQA2



Numéro de catalogue: 55072-1-AP

Informations de base

Numéro de catalogue:

55072-1-AP

Taille:

150ul, Concentration: 1000 µg/ml by Nanodrop and 507 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM_020056

Identification du gène (NCBI):

3118

Nom complet:

major histocompatibility complex, class II, DQ alpha 2

MW calculé

28 kDa

MW observés:

28 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:200-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, souris

Contrôles positifs:

WB : cellules Raji, tissu de thymus de souris, tissu splénique de souris

Informations générales

HLA-DQA2, also named as HLA-DXA, belongs to the MHC class II family. It binds peptides derived from antigens that access the endocytic route of antigen presenting cells (APC) and presents them on the cell surface for recognition by the CD4 T-cells. The antibody is specific to HLA-DQA2.

Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Raji cells were subjected to SDS PAGE followed by western blot with 55072-1-AP (HLA-DQA2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.