

À des fins de recherche uniquement

Anticorps Monoclonal anti-STRO-1

Numéro de catalogue: 65184-1-Ig



Informations de base

Numéro de catalogue:

65184-1-Ig

Taille:

100ug, 0.5 mg/ml

Hôte:

Mouse

Isotype:

IgM, lambda

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purifié par précipitation et chromatographie

CloneNo.:

STRO-1

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.



1X10⁶ lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.