

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD45RA

Numéro de catalogue: 65234-1-Ig



Informations de base

Numéro de catalogue:

65234-1-Ig

Taille:

100ug, 0.5 mg/ml

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

XM_006249910

Identification du gène (NCBI):

24699

Nom complet:

protein tyrosine phosphatase,
receptor type, C

MW calculé

143 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

OX-33

Applications

Applications testées:

FC

Spécificité de l'espèce:

rat

Informations générales

Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1×10^6 Wistar rat splenocytes were surface stained with 0.5 μ g Anti-Rat CD45RA (65234-1-Ig, Clone:OX-33) and FITC anti-Mouse IgG1 Antibody at dilution 1:100, or 0.5 μ g Control Antibody. Cells were not fixed.