

À des fins de recherche uniquement

Anticorps Monoclonal anti-FMNL1

Numéro de catalogue: 67587-1-Ig



Informations de base

Numéro de catalogue: 67587-1-Ig	Numéro d'acquisition GenBank: BC111564	Méthode de purification: Purification par protéine A
Taille: 150ul , Concentration: 1800 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 752	CloneNo.: 1B12A4
Hôte: Mouse	Nom complet: formin-like 1	Dilutions recommandées: WB 1:1000-1:6000
Isotype: IgG2b	MW observés: 140-150 kDa	
Immunogen Catalog Number: AG27322		

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules Jurkat, cellules HL-60, cellules K-562, tissu de thymus de porc
Spécificité de l'espèce: Humain, porc	

Informations générales

FMNL1, also named as C17orf1, C17orf1B and FMNL, belongs to the formin homology family. It may play a role in the control of cell motility and survival of macrophages. FMNL1 protein may be involved in anti-apoptosis mechanisms. It is overexpressed in lymphoid malignances and is found associated with Akt. 67587-1-Ig antibody detects the protein around 140-150 kDa which is similar to the results reported in publications. (PMID: 21868368, 22986496, 23182705)

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 67587-1-Ig (FMNL1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.

Jurkat cells were subjected to SDS PAGE followed by western blot with 67587-1-Ig (FMNL1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.

Various lysates were subjected to SDS PAGE followed by western blot with 67587-1-Ig (FMNL1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



NK-92 cells were subjected to SDS PAGE followed by western blot with 67587-1-Ig (FMNL1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.