

À des fins de recherche uniquement

Anticorps Monoclonal anti-REL

Numéro de catalogue: 67747-1-Ig



Informations de base

Numéro de catalogue: 67747-1-Ig	Numéro d'acquisition GenBank: NM_001291746	Méthode de purification: Purification par protéine G
Taille: 150ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 5966	CloneNo.: 1F6A4
Hôte: Mouse	Nom complet: v-rel reticuloendotheliosis viral oncogene homolog (avian)	Dilutions recommandées: WB 1:1000-1:6000
Isotype: IgG1	MW calculé 69 kDa	
Immunogen Catalog Number: AG26803	MW observés: 68 kDa, 80 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs:
Spécificité de l'espèce: Humain	WB : cellules Daudi, cellules HEK-293, cellules Jurkat, cellules MOLT-4

Informations générales

REL, also known as, Proto-oncogene c-Rel, is a 619 amino acid protein, which may play a role in differentiation and lymphopoiesis. NF-kappa-B is a homo- or heterodimeric complex formed by the Rel-like domain-containing proteins RELA/p65, RELB, NFkB1/p105, NFkB1/p50, REL and NFkB2/p52. The NF-kappa-B heterodimer RELA/p65-c-Rel is a transcriptional activator. The full-length 80-kDa c-Rel was identified (PMID: 9857058).

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 67747-1-Ig (REL antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.