

À des fins de recherche uniquement

Anticorps Monoclonal anti-Human IgG Kappa chain



Numéro de catalogue: 67761-1-Ig

Informations de base

Numéro de catalogue:

67761-1-Ig

Taille:

150ul, Concentration: 1000 µg/ml by Nanodrop and 439 µg/ml by Bradford method using BSA as the standard;

Hôte:

Mouse

Isotype:

IgG1

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

MW calculé

155 kDa

MW observés:

~25 kDa

Méthode de purification:

Purification par protéine A

CloneNo.:

4D11F7

Dilutions recommandées:

WB 1:5000-1:50000

IHC 1:2500-1:10000

Applications

Applications testées:

IHC, WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : plasma humain,

IHC : tissu d'amygdalite humain,

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9.0; (*) À défaut, 'le démasquage de l'antigène peut être 'effectué avec un tampon citrate pH 6.0.

Informations générales

This antibody detects the Kappa chain of human IgG,

Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:5000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:5000 (under 4x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:5000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



16 nL human plasma were subjected to SDS PAGE followed by western blot with 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.