

À des fins de recherche uniquement

Anticorps Monoclonal anti-TMEM74

Numéro de catalogue: 67773-1-Ig



Informations de base

Numéro de catalogue:	67773-1-Ig	Numéro d'acquisition GenBank:	BC030710	Méthode de purification:	Purification par protéine G
Taille:	150ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	157753	CloneNo.:	2H1B7
Hôte:	Mouse	Nom complet:	transmembrane protein 74	Dilutions recommandées:	WB 1:2000-1:10000
Isotype:	IgG1	MW calculé	305 aa, 33 kDa		
Immunogen Catalog Number:	AG22359	MW observés:	33 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : cellules LNCaP, cellules HEK-293, cellules L02
Spécificité de l'espèce:	Humain		

Informations générales

Transmembrane protein 74 (TMEM74), plays an essential role in autophagy during cell starvation and other stress conditions. TMEM74 is localized to the lysosome and autophagosome. It has been reported that TMEM74-induced autophagy may be associated with PI3K signal transduction. (PMID: 19029833, 18294959)

Stockage


Stockage:
Stocker à -20 °C.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



LNCaP cells were subjected to SDS PAGE followed by western blot with 67773-1-Ig (TMEM74 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.