

À des fins de recherche uniquement

Anticorps Monoclonal anti-PKMYT1

Numéro de catalogue: 67806-1-Ig **1 Publications**



Informations de base

Numéro de catalogue: 67806-1-Ig	Numéro d'acquisition GenBank: NML_004203	Méthode de purification: Purification par protéine A
Taille: 150ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 9088	CloneNo.: 1D12D3
Hôte: Mouse	Nom complet: protein kinase, membrane associated tyrosine/threonine 1	Dilutions recommandées: WB 1:5000-1:50000 IHC 1:500-1:2000
Isotype: IgG2a	MW calculé: 55 kDa	
Immunogen Catalog Number: AG30849	MW observés: 55 kDa	

Applications

Applications testées:

IHC, WB, ELISA

Demandes citées:

IHC, WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain

Contrôles positifs:

WB : cellules LNCaP, cellules 4T1, cellules HEK-293, cellules HeLa, cellules HSC-T6, cellules Jurkat, cellules K-562, cellules NIH/3T3, cellules PC-12

IHC : tissu de cancer du côlon humain,

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (*) À défaut, le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6,0.

Informations générales

Protein kinase, membrane associated tyrosine/threonine (PKMYT1) is a member of the Wee family of protein kinases. PKMYT1 acts as a negative regulator of the cell cycle by inactivating the CDK1-cyclinB complex through phosphorylation of Tyr14/Tyr15 and preventing cell cycle entering into mitosis at the G2/M transition. In mammalian somatic cells, PKMYT1 affects the assembly of the Golgi apparatus and endoplasmic reticulum during telophase of mitosis by inhibiting cyclin activities. Abnormal expression of PKMYT1 is closely associated with multiple types of cancers. PKMYT1 protein was mainly localized in the cytoplasm. (PMID: 31695486, PMID: 32801781)

Publications notables

Autrice	Pubmed ID	Journal	Application
Linhua Yang	36988891	Cell Oncol (Dordr)	IHC,WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 67806-1-Ig (PKMYT1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 67806-1-Ig (PKMYT1 antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).