

À des fins de recherche uniquement

Anticorps Monoclonal anti-SNX3

Numéro de catalogue: 68292-1-Ig



Informations de base

Numéro de catalogue:	68292-1-Ig	Numéro d'acquisition GenBank:	BC008444	Méthode de purification:	Purification par protéine G
Taille:	150ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	8724	CloneNo.:	2B5B7
Hôte:	Mouse	Nom complet:	sorting nexin 3	Dilutions recommandées:	WB 1:5000-1:50000 IF 1:400-1:1600
Isotype:	IgG1	MW calculé	19 kDa		
Immunogen Catalog Number:	AG1160	MW observés:	16-20 kDa		

Applications

Applications testées:	IF, WB, ELISA	Contrôles positifs:	
Spécificité de l'espèce:	Humain, rat, souris	WB :	cellules HeLa, cellules HepG2, cellules HSC-T6, cellules Jurkat, cellules K-562, cellules NIH/3T3
		IF :	cellules A431,

Informations générales

SNX3 is a member of the sorting nexin family of proteins regulating membrane traffic. SNX3 regulates endosomal function through its PX (phox homology)-domain-mediated interaction with PtdIns(3)P (PMID: 11433298). It is a potential regulator of neurite formation and may play an important role in the development and function of the central nervous system (PMID: 19576982; 20817026).

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 68292-1-Ig (SNX3 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using SNX3 antibody (68292-1-Ig, Clone: 2B5B7) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).