

À des fins de recherche uniquement

Anticorps Monoclonal anti-MED23

Numéro de catalogue: 68307-1-Ig



Informations de base

Numéro de catalogue:	68307-1-Ig	Numéro d'acquisition GenBank:	BC060759	Méthode de purification:	Purification par protéine G
Taille:	150ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	9439	CloneNo.:	1H4G4
Hôte:	Mouse	Nom complet:	mediator complex subunit 23	Dilutions recommandées:	WB 1:5000-1:50000
Isotype:	IgG1	MW calculé	156 kDa		
Immunogen Catalog Number:	AG28483	MW observés:	150 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : cellules A549, cellules 4T1, cellules HEK-293, cellules HeLa, cellules HSC-T6, cellules Jurkat, cellules K-562, cellules LNCaP, cellules NIH/3T3
Spécificité de l'espèce:	Humain, rat, souris		

Informations générales

MED23 is a component of the Mediator complex, it is a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional pre-initiation complex with RNA polymerase II and the general transcription factors.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 68307-1-Ig (MED23 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



4T1 cells were subjected to SDS PAGE followed by western blot with 68307-1-Ig (MED23 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.