

À des fins de recherche uniquement

Anticorps Recombinant de lapin anti-G3BP1



Numéro de catalogue: 80777-1-RR

Informations de base

Numéro de catalogue: 80777-1-RR	Numéro d'acquisition GenBank: BC006997	Méthode de purification: Purification par protéine A
Taille: 100ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 10146	CloneNo.: 7F14
Hôte: Lapin	Nom complet: GTPase activating protein (SH3 domain) binding protein 1	Dilutions recommandées: IF 1:500-1:2000
Isotype: IgG	MW calculé 466 aa, 52 kDa	
Immunogen Catalog Number: AG3728		

Applications

Applications testées: IF, ELISA	Contrôles positifs: IF : sodium arsenite treated HeLa cells,
Spécificité de l'espèce: Humain	

Informations générales

GAP SH3 Binding Protein 1 (G3BP1), also named as G3BP, is an effector of stress granule (SG) assembly. SG biology plays an important role in the pathophysiology of TDP-43 in ALS and FTLD-U. G3BP1 can be used as a marker of SG. It has been shown to function downstream of Ras and play a role in RNA metabolism, signal transduction, and proliferation. G3BP1 is a ubiquitously expressed protein that localizes to the cytoplasm in proliferating cells and to the nucleus in non-proliferating cells. G3BP1 has recently been implicated in cancer biology.

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed sodium arsenite treated HeLa cells using G3BP1 antibody (80777-1-RR, Clone: 7F14) at dilution of 1:1400 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).