

À des fins de recherche uniquement

Anticorps Monoclonal anti-FcεR1

Numéro de catalogue: **CL405-65091**



Informations de base

Numéro de catalogue:

CL405-65091

Taille:

100ug , 500 µg/ml

Hôte:

Armenian Hamster

Isotype:

IgG

Numéro d'acquisition GenBank:

BC125455

Identification du gène (NCBI):

14125

Nom complet:

Fc receptor, IgE, high affinity I, alpha polypeptide

Méthode de purification:

Purification par affinité

CloneNo.:

MAR-1

Excitation/Emission maxima

wavelengths:

399 nm / 422 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Fc epsilon Receptor I alpha (FcεRIα) is a single-pass type I membrane protein which contains 2 immunoglobulin-like domains. FcεRIα forms a tetrameric complex with one beta and two gamma subunits (PMID: 19406478). FcεRI is the high-affinity receptor for the Fc region of immunoglobulin E (IgE) (PMID: 1532373). It is expressed on mast cells and basophils and plays a central role in the allergic response (PMID: 2148316).

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ MC/9 cells were surface stained with 1 ug CoraLite® Plus 405 Anti-Mouse FcεR1 (CL405-65091, Clone:MAR-1) (red), or 1 ug Isotype Control. Cells were not fixed.