

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD94

Numéro de catalogue: **CL488-65240**



Informations de base

Numéro de catalogue:
CL488-65240

Taille:
100ug , 0.5 mg/ml

Hôte:
Rat

Isotype:
IgG2a, kappa

Numéro d'acquisition GenBank:
BC120854

Identification du gène (NCBI):
16643

Nom complet:
**killer cell lectin-like receptor,
subfamily D, member 1**

Méthode de purification:
N/A

CloneNo.:
18D3

Excitation/Emission maxima
wavelengths:
493 nm / 522 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
souris

Informations générales


Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 488 Anti-Mouse CD94 (CL488-65240, Clone:18D3), or 0.5 ug Isotype Control. Cells were not fixed.