

À des fins de recherche uniquement

Anticorps Monoclonal anti-PDCD4

Numéro de catalogue: **CL488-66100** **Phare**



Informations de base

Numéro de catalogue:	Numéro d'acquisition GenBank:	Méthode de purification:
CL488-66100	BC026104	Purification par protéine G
Taille:	Identification du gène (NCBI):	CloneNo.:
100ul , Concentration: 1000 µg/ml by Nanodrop;	27250	1F1A11
Hôte:	Nom complet:	Dilutions recommandées:
Mouse	programmed cell death 4 (neoplastic transformation inhibitor)	WB 1:500-1:1000
Isotype:	MW calculé	Excitation/Emission maxima wavelengths:
IgG1	469 aa, 52 kDa	493 nm / 522 nm
Immunogen Catalog Number:	MW observés:	
AG19239	56 kDa	

Applications

Applications testées:	Contrôles positifs:
WB	WB : cellules Jurkat, cellules HEK-293, cellules HeLa
Spécificité de l'espèce:	
Humain	

Informations générales

Programmed cell death 4 (Pdc4) is a novel tumor suppressor that inhibits translation, progression and invasion. It was first identified as being differentially upregulated during apoptosis. Pdc4 interferes with the activity of the eukaryotic initiation factor (eIF) 4A by displacing the scaffold protein eIF4G from its binding to the RNA helicase eIF4A.

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with CL488-66100 (PDCD4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.