

À des fins de recherche uniquement

Anticorps Monoclonal anti-SNX3

Numéro de catalogue: **CL488-68292**



Informations de base

Numéro de catalogue:	Numéro d'acquisition GenBank:	Méthode de purification:
CL488-68292	BC008444	Purification par protéine G
Taille:	Identification du gène (NCBI):	CloneNo.:
100ul , Concentration: 1000 µg/ml by Nanodrop;	8724	2B5B7
Hôte:	Nom complet:	Dilutions recommandées:
Mouse	sorting nexin 3	IF 1:50-1:500
Isotype:	MW calculé	Excitation/Emission maxima wavelengths:
IgG1	19 kDa	493 nm / 522 nm
Immunogen Catalog Number:	MW observés:	
AG1160	16-20 kDa	

Applications

Applications testées:	Contrôles positifs:
IF	IF : cellules A431,
Spécificité de l'espèce:	
Humain, rat, souris	

Informations générales

SNX3 is a member of the sorting nexin family of proteins regulating membrane traffic. SNX3 regulates endosomal function through its PX (phox homology)-domain-mediated interaction with PtdIns(3)P (PMID: 11433298). It is a potential regulator of neurite formation and may play an important role in the development and function of the central nervous system (PMID: 19576982; 20817026).

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using CoraLite® Plus 488 SNX3 antibody (CL488-68292, Clone: 2B5B7) at dilution of 1:200.