

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD8a

Numéro de catalogue: **CL555-65205**



Informations de base

Numéro de catalogue:
CL555-65205

Taille:
100ug, 0.5 mg/ml

Hôte:
Rat

Isotype:
IgG2b, lambda

Numéro d'acquisition GenBank:
BC030679

Identification du gène (NCBI):
12525

Nom complet:
CD8 antigen, alpha chain

Méthode de purification:
Purification par affinité

CloneNo.:
5H10-1

Excitation/Emission maxima
wavelengths:
554 nm / 570 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
souris

Informations générales

CD8 is a transmembrane glycoprotein composed of two disulfide-linked chains. It can be present as a homodimer of CD8a or as a heterodimer of CD8a and CD8β (PMID: 3264320; 8253791). CD8 is found on most thymocytes. The majority of class I-restricted T cells express mostly the CD8aβ heterodimer while CD8aa homodimers alone have been found on some gut intraepithelial T cells, on some T cell receptor (TCR) γδ T cells and on NK cells (PMID: 2111591; 1831127; 8420975). CD8 acts as a co-receptor that binds to MHC class-I and participates in cytotoxic T cell activation (PMID: 8499079). During T cell development, CD8 is required for positive selection of CD4-/CD8+ T cells (PMID: 1968084).

Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface co-stained with APC Anti-Mouse CD3 and 0.5 ug CoraLite® Plus 555 Anti-Mouse CD8a (CL555-65205, Clone:5H10-1) or 0.5 ug Isotype Control. Cells were not fixed.