

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD49b

Numéro de catalogue: **CL555-65252**



Informations de base

Numéro de catalogue:
CL555-65252

Taille:
100 tests, 5 µl/test

Hôte:
Mouse

Isotype:
IgG1, kappa

Numéro d'acquisition GenBank:
NM_002203

Identification du gène (NCBI):
3673

Nom complet:
integrin, alpha 2 (CD49B, alpha 2
subunit of VLA-2 receptor)

MW calculé
126 kDa

Méthode de purification:
Purification par affinité

CloneNo.:
AK7

Excitation/Emission maxima
wavelengths:
554 nm / 570 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
Humain

Informations générales

Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 555 Anti-Human CD49b (CL555-65252, Clone: AK7) (red) or unstained. Cells were not fixed. Monocytes were gated.