

À des fins de recherche uniquement

# Anticorps Monoclonal anti-Mouse IgG2b isotype control



Numéro de catalogue: CL555-66360-3

## Informations de base

Numéro de catalogue:

CL555-66360-3

Taille:

100ug, Concentration: 0.2 mg/mL by Nanodrop;

Hôte:

Mouse

Isotype:

IgG2b

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par protéine A

CloneNo.:

11B8C4

Excitation/Emission maxima wavelengths:

557 nm / 570 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

Megathura crenulata

## Informations générales

This antibody is used as the isotype control of mouse IgG2b. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

## Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug CoraLite®555 Anti-Mouse IgG2b isotype control (CL555-66360-3, Clone:11B8C4) (red), or CoraLite®555 Anti-Human CD4 (OKT4) (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.