

À des fins de recherche uniquement

Anticorps Monoclonal anti-SGTA

Numéro de catalogue: **CL594-60305**



Informations de base

Numéro de catalogue: CL594-60305	Numéro d'acquisition GenBank: BC008885	Méthode de purification: Purification par protéine A
Taille: 100ul , Concentration: 2000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 6449	CloneNo.: 3B10B11
Hôte: Mouse	Nom complet: small glutamine-rich tetrapeptide repeat (TPR)- containing, alpha	Dilutions recommandées: IF 1:50-1:500
Isotype: IgG2b	MW calculé 34 kDa	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: AG1475	MW observés: 34-37 kDa	

Applications

Applications testées: IF	Contrôles positifs: IF : cellules HeLa,
Spécificité de l'espèce: Humain, souris	

Informations générales

Small glutamine-rich tetrapeptide repeat-containing protein alpha (SGTA) is a steroid receptor molecular co-chaperone critical in a range of biological processes, including viral infection, cell division, mitosis, and cell cycle checkpoint activation. Overexpression of SGTA has been found in various cancers, including prostate/ovary/testis/liver/lung cancer.

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CL594-60305 (SGTA antibody) at dilution of 1:100. .