

À des fins de recherche uniquement

Anticorps Monoclonal anti-DNAJB6

Numéro de catalogue: CL594-66587



Informations de base

Numéro de catalogue: CL594-66587	Numéro d'acquisition GenBank: BC000177	Méthode de purification: Purification par protéine G
Taille: 100ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 10049	CloneNo.: 1A1D2
Hôte: Mouse	Nom complet: DnaJ (Hsp40) homolog, subfamily B, member 6	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Isotype: IgG1	MW calculé 241 aa, 27 kDa	
Immunogen Catalog Number: AG2258	MW observés: 27 kDa, 38 kDa	

Applications

Applications testées:
FC (Intra)

Spécificité de l'espèce:
Humain, porc, rat, souris

Informations générales

DNAJB6, belongs to the heat shock protein 40 (HSP40) family, a noncanonical member of the DNAJ-chaperone family, is mainly expressed in neurons. It plays an important role in recovery from misfolded protein aggregates (PMID: 27276715). As seen by the western blot stained for DNAJB6, there are 2 splice forms of the protein resulting in a longer isoform (38 kDa) and a shorter one (27 kDa). DNAJB6 (27 kDa) isoform but not DNAJB6 (38 kDa) isoform is effective with suppressing cytoplasmic protein aggregation (PMID: 20159555). The antibody 66587-1-Ig recognizes both isoforms.

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.

Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ Jurkat cells were intracellularly stained with 0.4 ug CoraLite® 594 Anti-Human DNAB6 (CL594-66587, Clone:1A1D2) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).