

À des fins de recherche uniquement

Anticorps Monoclonal anti-BCAR1

Numéro de catalogue: **CL594-67215**



Informations de base

Numéro de catalogue:	CL594-67215	Numéro d'acquisition GenBank:	BC062556	Méthode de purification:	Purification par protéine A
Taille:	100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	9564	CloneNo.:	2D11C10
Hôte:	Mouse	Nom complet:	breast cancer anti-estrogen resistance	Dilutions recommandées:	IF 1:50-1:500
Isotype:	IgG2a	MW calculé	870 aa, 93 kDa	Excitation/Emission maxima wavelengths:	588 nm / 604 nm
Immunogen Catalog Number:	AG10604	MW observés:	116 kDa		

Applications

Applications testées:	FC (Intra), IF	Contrôles positifs:	IF : tissu de cancer du sein humain,
Spécificité de l'espèce:	Humain		

Informations générales

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed human breast cancer tissue using CoraLite®594 BCAR1 antibody (CL594-67215, Clone: 2D11C10) at dilution of 1:100.



Immunofluorescent analysis of (4% PFA) fixed human breast cancer tissue using CoraLite®594 BCAR1 antibody (CL594-67215, Clone: 2D11C10) at dilution of 1:100.



1X10⁶ A431 cells were intracellularly stained with 0.4 ug CoraLite®594 Anti-Human BCAR1 (CL594-67215, Clone:2D11C10) (red), or 0.4 ug Mouse IgG2a Isotype Control (CL594-66360-2, Clone: K11A1B2A2) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).