

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD41

Numéro de catalogue: **CL647-65189**



## Informations de base

Numéro de catalogue:

CL647-65189

Taille:

100ug, 500 µg/ml

Hôte:

Rat

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC120493

Identification du gène (NCBI):

16399

Nom complet:

integrin alpha 2b

Méthode de purification:

Purification par affinité

CloneNo.:

MWRReg30

Excitation/Emission maxima

wavelengths:

654 nm / 674 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

## Informations générales

### Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> Balb/c mouse peripheral blood platelets were surface stained with 0.5 ug Coralite® Plus 647 Anti-Mouse CD41 (CL647-65189, Clone:MWRReg30), or 0.5 ug Isotype Control. Cells were not fixed.



1X10<sup>6</sup> Balb/c mouse peripheral blood platelets were surface stained with 0.5 ug Coralite® Plus 647 Anti-Mouse CD41 (CL647-65189, Clone:MWRReg30) (red), or 0.5 ug Isotype Control. Cells were not fixed.