

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD48

Numéro de catalogue: **CL647-65227**



## Informations de base

Numéro de catalogue:  
CL647-65227

Taille:  
100ug, 0.5 mg/ml

Hôte:  
Armenian Hamster

Isotype:  
IgG

Numéro d'acquisition GenBank:  
BC060977

Identification du gène (NCBI):  
12506

Nom complet:  
CD48 antigen

Méthode de purification:  
Purification par protéine A

CloneNo.:  
HM48-1

Excitation/Emission maxima  
wavelengths:  
654 nm / 674 nm

## Applications

Applications testées:  
FC

Spécificité de l'espèce:  
souris

## Informations générales

### Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Données de validation sélectionnées



1X10<sup>6</sup> mouse splenocytes were surface stained with 1 ug CoraLite® Plus 647 Anti-Mouse CD48 (CL647-65227, Clone: HM48-1) (red), or 1 ug Control Antibody. Cells were not fixed.