

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD45RA

Numéro de catalogue: **CL647-65234**



Informations de base

Numéro de catalogue:
CL647-65234

Taille:
100ug , 0.5 mg/ml

Hôte:
Mouse

Isotype:
IgG1, kappa

Numéro d'acquisition GenBank:
XM_006249910

Identification du gène (NCBI):
24699

Nom complet:
protein tyrosine phosphatase,
receptor type, C

MW calculé
143 kDa

Méthode de purification:
Purification par affinité

CloneNo.:
OX-33

Excitation/Emission maxima
wavelengths:
654 nm / 674 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
rat

Informations générales


Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 647 Anti-Rat CD45RA (CL647-65234, Clone:OX-33), or 0.5 ug Control Antibody. Cells were not fixed.