

À des fins de recherche uniquement

Anticorps Monoclonal anti-WASF3

Numéro de catalogue: **CL647-67620**



Informations de base

Numéro de catalogue: CL647-67620	Numéro d'acquisition GenBank: BC050283	Méthode de purification: Purification par protéine G
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 10810	CloneNo.: 3B7D10
Hôte: Mouse	Nom complet: WAS protein family, member 3	Dilutions recommandées: IF 1:50-1:500
Isotype: IgG1	MW calculé 55 kDa	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Immunogen Catalog Number: AG29991	MW observés: 50-55 and 66-70 kDa	

Applications

Applications testées: IF	Contrôles positifs: IF : tissu cérébral de souris,
Spécificité de l'espèce: Humain, Lapin, porc, rat, souris	

Informations générales

WASF3, also named as WAVE3 or SCAR3, belongs to the SCAR/WAVE family. It plays a role in the regulation of cell morphology and cytoskeletal organization. WASF3 is associated with the movement, invasion and metastasis in different cancers such as breast cancer, gastric cancer and prostate cancer. It has been shown that WASF3 levels are elevated in advanced stage cancers compared to lower grade tumors and normal tissue. 67620-1-Ig antibody recognizes 50-55 kDa and 66-70 kDa protein similar to the papers.(PMID: 21544801, 15826941, 31038356, 27432794, 31542393)

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using CoraLite® Plus 647 WASF3 antibody (CL647-67620, Clone: 3B7D10) at dilution of 1:200.