

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD41

Numéro de catalogue: **CL750-65189**



Informations de base

Numéro de catalogue:

CL750-65189

Taille:

100ug, 0.5 mg/ml

Hôte:

Rat

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC120493

Identification du gène (NCBI):

16399

Nom complet:

integrin alpha 2b

Méthode de purification:

Purification par affinité

CloneNo.:

MWRReg30

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ Balb/c mouse peripheral blood platelets were surface stained with 0.2 ug CoraLite® Plus 750 Anti-Mouse CD41 (CL750-65189, Clone: MWRReg30) (red), or 0.2 ug Isotype Control. Cells were not fixed.



1X10⁶ Balb/c mouse peripheral blood platelets were surface stained with 0.2 ug CoraLite® Plus 750 Anti-Mouse CD41 (CL750-65189, Clone: MWRReg30), or 0.2 ug Isotype Control. Cells were not fixed.