

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD11b

Numéro de catalogue: **CL750-65229**



Informations de base

Numéro de catalogue:

CL750-65229

Taille:

100ug, 0.5 mg/ml

Hôte:

Mouse

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

NM_012711

Identification du gène (NCBI):

25021

Nom complet:

integrin alpha M

MW calculé

127 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

OX42

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

rat

Informations générales

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). CD11b, also known as Integrin alpha M or CR3A, belongs to the integrin alpha chain family. CD11b forms an α/β heterodimer with CD18 (integrin $\beta 2$). CD11b/CD18 is implicated in various adhesive interactions of monocytes, macrophages and granulocytes as well as in mediating the uptake of complement-coated particles and pathogens (PMID: 9558116; 20008295). CD11b/CD18 is a receptor for the complement protein fragment iC3b, and is also a receptor for fibrinogen, factor X and ICAM1 (PMID: 2971974; 15485828).

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Rat CD11b (CL750-65229, Clone: OX42) or 0.5 ug Control Antibody. Cells were not fixed.